

# LABORHITS

## Midi Centrifuge G-L

For all types of separation processes

High-performance motor



### Features

- With metal cover
- Starts only when the cover is firmly closed
- Direct and easy setting of speed, rotational speed and programmes via membrane keypad

Due to the material separation at medium rotational speeds, this midi centrifuge is ideal for all types of separation processes, such as in molecular biology. It can be used with standard 1.5 / 2 ml tubes as well as with PCR strips. The cover opens automatically at the press of a button.

**Scope of delivery:** Midi centrifuge incl. metal rotor and cover for 12 vessels, with a capacity of 1.5/2 ml

### Specifications

**Max. RCF:** 16,500 x g  
**Max. speed:** 15,700 min<sup>-1</sup>  
**Max. capacity:** 12 x 1.5 / 2.0 ml  
**Power consumption:** 65 W  
**Dimensions (W x D x H):** 200 x 280 x 140 mm  
**Weight:** 3.9 kg  
**Power supply:** 220 ... 230 V, 50 ... 60 Hz  
**Protection class:** IP 20

Type	PK	Cat. No.	Price €	Special price €
G-L	1	4.664 755	999,00	699,00

All prices exclusive VAT or other local taxes.  
 Special price valid until 31.03.2020.